

PLACE



PLATFORM OF LABORATORIES FOR ADVANCES IN CARDIAC EXPERIENCE

ROMA 2024

- 12 Giugno **MINICORSI** Precongressuali
- 13-15 Giugno **CONGRESSO**

Prevenzione dello Stroke Cardioembolico

Live in the box:

Chiusura percutanea di auricolare sinistra

Dr. Enrico Mura

Direttore Cardiologia Interventistica HSF



ASL 3 Nuoro




Contrastless-no touch LAAO



Donna di 83 anni

- 04/2024: Stroke ischemico sotto terapia anticoagulante orale (Apixaban 2,5 mg x 2/die)
- IRC: filtrato glomerulare 19 ml/min/mq
- Fibrillazione atriale permanente
- Anemizzazioni ricorrenti (Angiodisplasie intestinali)












Table 1: the device sizing chart for single lobe LAA

LZ(mm)	LAmbre™ Occluder	[D] Umbrella Diameter(mm)	[D1] Cover Diameter(mm)	LAmbre™ Delivery System
11.0~12.5	LT-LAA-1622	16	22	8/9/10Fr; 900mm
12.5~14.0	LT-LAA-1824	18	24	
14.0~16.0	LT-LAA-2026	20	26	9/10Fr; 900mm
16.0~18.0	LT-LAA-2228	22	28	
18.0~20.0	LT-LAA-2430	24	30	
20.0~22.0	LT-LAA-2632	26	32	10Fr; 900mm
22.0~24.0	LT-LAA-2834	28	34	
24.0~26.0	LT-LAA-3036	30	36	
24.0~26.0	LT-LAA-3236	32	36	
26.0~28.0	LT-LAA-3438	34	38	
28.0~30.0	LT-LAA-3640	36	40	





Table 2: the device sizing chart for double or multiple lobes LAA



LZ(mm)	LAmbre™ Occluder	[D] Umbrella Diameter(mm)	[D1] Cover Diameter(mm)	LAmbre™ Delivery System
11.0~12.5	LT-LAA-1630	16	30	9/10Fr; 900mm
12.5~14.0	LT-LAA-1832	18	32	
14.0~16.0	LT-LAA-2032	20	32	10Fr; 900mm
16.0~18.0	LT-LAA-2234	22	34	
18.0~20.0	LT-LAA-2436	24	36	
20.0~22.0	LT-LAA-2638	26	38	


Note:

1. LZ: Landing zone
2. Please make sure the depth from the LZ to the distal end of the LAA ≥ 10 mm
3. Please make sure the Umbrella is deployed beyond the left circumflex artery














Grazie per l'attenzione